

TEAR

For Emergency Temporary External Use Only

**CELOX<sup>TM</sup> RAPID****Hemostatic Z-Fold Gauze**Chito-R<sup>TM</sup> Technology, Stops Bleeding Fast\*

To Control Moderate-To-Severe Bleeding

Works in Hypothermic Conditions\*

One Gauze Strip 3in x 5ft (7.6cm x 1.5m)

Do not  
re-useDo not  
resterilize**STERILE R**Sterilized using  
irradiation

Caution

Do not use if  
package is  
damaged**Rx only****Instructions for Use on Back**

\*Proven in laboratory trials to stop major arterial bleeding by coagulating the blood within 60 seconds, reducing blood loss. Granules shown to be effective in blood temperatures down to 56°F (hypothermic).  
Data on file at Medtrade Products Ltd.

TEAR

TEAR

TEAR

**FOR USE BY TRAINED EMERGENCY RESPONDERS****Instructions for Use:**

**1.** Identify the point of bleeding and apply direct pressure. Swab excess blood where practical.



**2.** Tear open pouch. Take out the CELOX RAPID and take one end of Z-folded gauze.



**3.** Tightly pack the unfolding CELOX RAPID directly to the source of bleeding, feeding in the CELOX RAPID until the wound is filled above skin level. Excess CELOX RAPID can be cut or torn if necessary.



**4.** Apply FIRM pressure directly to wound for 60 seconds. If bleeding persists apply additional pressure until bleeding stops.



**5.** Secure with a dressing to maintain pressure on the wound. Discard any remaining CELOX RAPID.

**6.** Transfer patient to Medical Facility as soon as possible.

**7.** Show empty pack to medical personnel.

**ATTENTION MEDICAL FACILITY PERSONNEL**

1. Physically remove CELOX RAPID from the wound and any loose surface granules.
2. Fully flood entire wound area with sterile saline irrigation solution.
3. Proceed with normal irrigation and/or suction.
4. Ensure all product is removed from the wound prior to initiation of wound treatment.

**Warnings:** Device is not intended for surgical use. Do not use in abdominal wounds and wounds not amenable to pressure. Loss of sterility potentially poses a risk of infection. Do not apply over eyes. If eye irritation occurs flush with water for 5 minutes. Contains chitosan from shellfish. Allergy studies show no adverse reaction to date. Data on file at Medtrade Products Ltd.

Manufactured by:  
Medtrade Products Ltd, Electra House,  
Crewe Business Park, Crewe, CW1 6GL, UK.  
[www.celoxmedical.com](http://www.celoxmedical.com)

Issue Date: 20 Feb 2020. Artwork Ref: MTP-19-1422

**REF** FG08839011

TEAR

TEAR